

Stabilization Using PEEK Rods - Clinical Outcomes in Patients with Grade 1 Degenerative Lumbar Spondylolisthesis

AUTHORS: Dr. Pablo Bouillot, Hôpital privé les Franciscaines, Nimes, France*
Sylvie Poulette, PhD, CEO at conseil et gestion de projets (clinique & santé) (Consulting and Project Management (Clinic & Health)) Bourdeaux, France**

ABSTRACT: Originates from Stabilization Using Peek Rods - Clinical outcomes in patients with grade 1 degenerative lumbar spondylolisthesis White Paper, Invibio Biomaterial Solutions™, 2018

LEVEL OF EVIDENCE: Level 2a Retrospective Cohort Study

Introduction

Instrumented rigid arthrodesis is commonly used for the surgical treatment of lumbar degenerative diseases. Nevertheless, rigid materials may contribute to stress shielding in the anterior column leading to risk of implant failure and adjacent segment disease (ASD),¹⁻¹⁰ responsible for a high rate of further surgery. Polyether ether ketone (PEEK) provides a modulus of elasticity similar to that of bone¹¹ that may limit these risks. Comparative biomechanical studies have shown that PEEK systems provide intervertebral stability comparable to rigid constructs, while allowing for a better redistribution of segmental loads and a reduced stress at the bone-screw interface.¹²⁻¹⁵ This retrospective study was carried out to evaluate the efficacy and the safety of the flexible stabilization system Initial VEOS PEEK®-Optima (Innov'Spine, France) without arthrodesis in patients treated for grade 1 degenerative lumbar spondylolisthesis.

Methods

Out of a homogeneous cohort of 66 patients operated with the flexible stabilization Initial VEOS PEEK-Optima Osteosynthesis System (made with PEEK-OPTIMA™ polymer from Invibio Biomaterial Solutions, Figure 1) without

spinal pathology, vertebral fracture, severe osteoporosis, neoplasia, inflammatory rheumatologic disease or active psychiatric illness were not included). Surgery was performed using PEEK-OPTIMA™ Rods and Pedicle Screws (Initial VEOS PEEK-Optima Osteosynthesis System), without arthrodesis and was associated with a foraminotomy. The patients' follow-up ranged from 1 to 4 years (mean 24 months). Clinical and radiological outcomes were assessed using the Oswestry Disability Index (ODI), complications, CT-scan, dynamic and static X-rays. Patient satisfaction was also investigated using the Patient Satisfaction Index (PSI).

Results

The efficacy of the treatment was evaluated by comparing the preoperative and postoperative ODI. The average ODI decreased from 50.53% to 8.00% and the rate of patients with moderate to total disability (ODI >15) from 89.5% (N=34) to 5.3% (N=2) after 1 to 4 years of follow-up. At the end of follow-up, 60.5% of the patients (N=23 patients) had no more disability and no patient (0%) had severe or total disability (compared to 34.2%, N=13 before surgery). Results are summarized in Table 1 and Figure 2 illustrates the improvement of pain-disability after surgery. Static

Treatment Efficacy: ODI Value Comparison

ODI Range	Before Surgery N = 38 (%)	After Surgery N = 38 (%)
0 – 4: no disability	0 (0%)	23 (60.49%)
5 – 14: minimal disability	4 (10.52%)	13 (34.19%)
15 – 24: moderate disability	21 (55.23%)	2 (5.26%)
25 – 34: severe disability	10 (26.3%)	0 (0%)
> 34: total disability	3 (7.89%)	0 (0%)

Table 1: Distribution in degrees of severity of pain-disability before and after surgery (1 to 4 years following surgery).

and dynamic radiographs showed a complete integration of the material without breaking or mobilization. Mobility on dynamic flexion/extension X-rays was observed in 63.2% of cases (N=24) and all the patients showed a perfect tolerance to the flexion/extension mobilization on clinical examination. Two (2) complications (one superficial sepsis and one misplaced screw) were associated with the



Image provided courtesy of Innov'Spine

Figure 1: The flexible stabilization system INITIAL VEOS PEEK®-Optima Osteosynthesis System (A) with PEEK-OPTIMA Rods and Pedicle Screws (B). PEEK rods provide elastic modulus similar to that of bone that offers both adequate rigidity for fusion and a flexibility that limits the stress created by rigid rods (C).

arthrodesis for treatment of grade 1 degenerative lumbar spondylolisthesis associated with severe canal stenosis, 38 were included in the study (10 patients lost to follow-up and 18 patients with neurological disease, other




Figure 2: Improvement of the pain-disability after surgery (1 to 4 years following surgery).

surgical procedure and required an additional surgical procedure. Clinical outcomes were supported by patients' satisfaction, using the PSI questionnaire, with a majority of patients (84%) quite satisfied with the treatment (Figure 3).




Figure 3: Patient satisfaction evaluated using the PSI questionnaire (answers to questionnaires collected between 1 to 4 years of follow-up)

Conclusion

The flexible Initial VEOS PEEK-Optima Osteosynthesis System, made from PEEK-OPTIMA Natural polymer, used as a rods-screws construct alone (without arthrodesis or instrumented interbody fusion) provides a safe and effective stabilization showing substantial improvement in ODI and a low rate of serious complications and reoperations, that may be considered as an alternative solution to rigid fusion systems in the treatment of grade 1 degenerative lumbar spondylolisthesis. Long-term, multicenter and comparative studies are nevertheless needed to demonstrate improved benefits over rigid fixation. ▲

ABOUT THE AUTHORS

Dr. Pablo Bouillot

Dr. Pablo Bouillot, Neurosurgeon, Hôpital privé les Franciscaines, Nîmes, France.



Sylvie Poulette, PhD

Sylvie Poulette, PhD, CEO at conseil et gestion de projets (clinique & santé), (Consulting and Project Management (Clinic & Health)) Bourdeaux, France. Dr. Poulette holds a PhD in Biological and Medical Sciences.



REFERENCES

- Lehmann, T. R. et al. Long-term follow-up of lower lumbar fusion patients. *Spine (Phila. Pa. 1976)*. (1987).
- Ghiselli, G., Wang, J. C., Bhatia, N. N., Hsu, W. K. & Dawson, E. G. Adjacent segment degeneration in the lumbar spine. *J. Bone Jt. Surg. - Ser. A* (2004).
- Ren, C., Song, Y., Liu, L. & Xue, Y. Adjacent segment degeneration and disease after lumbar fusion compared with motion-preserving procedures: a meta-analysis. *Eur. J. Orthop. Surg. Traumatol.* (2014).
- Rahm, M. D. & Hall, B. B. Adjacent-segment degeneration after lumbar fusion with instrumentation: a retrospective study. *Journal of spinal disorders* (1996).
- Lee, C. K. Accelerated degeneration of the segment adjacent to a lumbar fusion. *Spine (Phila. Pa. 1976)*. (1988).
- Eck, J. C., Humphreys, S. C. & Hodges, S. D. Adjacent-segment degeneration after lumbar fusion: a review of clinical, biomechanical, and radiologic studies. *Am. J. Orthop. (Belle Mead, NJ)*. (1999).
- Hilibrand, A. S. & Robbins, M. Adjacent segment degeneration and adjacent segment disease: The consequences of spinal fusion? *Spine J.* (2004).
- Park, P., Garton, H. J., Gala, V. C., Hoff, J. T. & McGillicuddy, J. E. Adjacent segment disease after lumbar or lumbosacral fusion: Review of the literature. *Spine* (2004).
- Harrop, J. S. et al. Lumbar adjacent segment degeneration and disease after arthrodesis and total disc arthroplasty. *Spine (Phila. Pa. 1976)*. (2008).
- Lee, C. S. et al. Risk factors for adjacent segment disease after lumbar fusion. *Eur. Spine J.* (2009).

11. Vadapalli, S. et al. Biomechanical rationale for using polyetheretherketone (PEEK) spacers for lumbar interbody fusion - A finite element study. *Spine (Phila. Pa. 1976)*. (2006).
12. Turner, J. L., Paller, D. J. & Murrell, C. B. The mechanical effect of commercially pure titanium and polyetheretherketone rods on spinal implants at the operative and adjacent levels. *Spine (Phila. Pa. 1976)*. (2010).
13. Ponnappan, R. K. et al. Biomechanical evaluation and comparison of polyetheretherketone rod system to traditional titanium rod fixation. *Spine Journal* (2009).
14. Bruner, H. J. et al. Biomechanics of polyaryletherketone rod composites and titanium rods for posterior lumbosacral instrumentation. *J. Neurosurg. Spine* (2010).
15. Chou, W. K., Chien, A. & Wang, J. L. Biomechanical analysis between PEEK and titanium screw-rods spinal construct subjected to fatigue loading. *J. Spinal Disord. Tech.* (2015).

* Since 2017, Dr. Pablo Bouillot, has provided ad hoc consultancy services to Invibio Ltd.

* The abstract and conclusions presented have been provided by a practicing orthopedic surgeon. His view and experiences are his own and do not necessarily reflect those of others. "Invibio" disclaims any liabilities or loss in connection with the information herein.

** Since 2017, Sylvie Poulette, PhD has provided ad hoc consultancy services to Invibio Ltd.

The third party trademarks used herein are the trademarks of their respective owners.

Copyright ©2018 Invibio Ltd. INVIBIO™, PEEK-OPTIMA™, INVIBIO BIOMATERIAL SOLUTIONS™ are trademarks of Victrex plc or its group companies. All rights reserved.



Invibio Ltd.

Victrex Technology Centre
Hillhouse International
Thornton-Cleveleys
Lancashire
FY5 4QD, UK

Tel: +44 (0) 1253 898 000

FAX: +44 (0) 1253 898 001

Invibio Inc.

300 Conshohocken State Road
West Conshohocken, PA
19428
USA

Toll Free: 866-INVIBIO (468-4246)

Tel: (484) 342-6004

Fax: (484) 342-6005

**For further information please email us at info@invibio.com
or visit our website at:**

► Invibio.com

Victrex plc and/or its group companies ("Victrex plc") believes that the information in this document is an accurate description of the typical characteristics and/or uses of the product or products, but it is the customer's responsibility to thoroughly test the product in each specific application to determine its performance, efficacy, and safety for each end-use product, device or other application. Suggestions of uses should not be taken as inducements to infringe any particular patent. The information and data contained herein are based on information we believe reliable. Mention of a product in this document is not a guarantee of availability.

Victrex plc reserves the right to modify products, specifications and/or packaging as part of a continuous program of product development. Victrex plc makes no warranties, express or implied, including, without limitation, a warranty of fitness for a particular purpose or of intellectual property non-infringement, including, but not limited to patent non-infringement, which are expressly disclaimed, whether express or implied, in fact or by law.

Further, Victrex plc makes no warranty to your customers or agents, and has not authorized anyone to make any representation or warranty other than as provided above. Victrex plc shall in no event be liable for any general, indirect, special, consequential, punitive, incidental or similar damages, including without limitation, damages for harm to business, lost profits or lost savings, even if Victrex has been advised of the possibility of such damages regardless of the form of action.

Supporting information is available on request for all claims referenced in this document.

Copyright ©2018 Invibio Ltd. INVIBIO™, JUVORA™ PEEK-OPTIMA™, INVIBIO BIOMATERIAL SOLUTIONS™ are trademarks of Victrex plc or its group companies. All rights reserved.